

Product datasheet for **TP725086**

CYP26C1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytochrome P450 26C1
Species:	Human
Expression cDNA Clone or AA Sequence:	F230-L360
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cytochrome P450 26C1 is produced by E.coli expression system and the target gene encoding F230-L360 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	340665
UniProt ID:	Q6V0L0



[View online »](#)