

Product datasheet for **TP725089**

TUBB2B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tubulin beta-2B chain
Species:	Human
Expression cDNA Clone or AA Sequence:	C12-E205
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Tubulin beta-2B chain is produced by E.coli expression system and the target gene encoding C12-E205 is expressed.
Stability:	12 months from date of despatch
Locus ID:	347733
UniProt ID:	Q9BVA1



[View online »](#)