

Product datasheet for **TP725090**

Myosin 1G (MYO1G) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Unconventional myosin-Ig/MYO1G
Species:	Human
Expression cDNA Clone or AA Sequence:	R190-I391
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Unconventional myosin-Ig/MYO1G is produced by E.coli expression system and the target gene encoding R190-I391 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	64005
UniProt ID:	B0I1T2
Synonyms:	Minor histocompatibility antigen HA-2; mHag HA-2; HA2



[View online »](#)