

Product datasheet for **TP725092**

EEF1E1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human EEF1E1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-H174
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human EEF1E1 Protein is produced by E.coli expression system and the target gene encoding M1-H174 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9521
UniProt ID:	O43324
Synonyms:	Eukaryotic translation elongation factor 1 epsilon-1; Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3; Elongation factor p18; Multisynthase complex auxiliary component p18 ; AIMP3; P18
Protein Families:	Druggable Genome



[View online »](#)