

Product datasheet for **TP725098**

HNRPH3 (HNRNPH3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human hnRNP H3(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Y346
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human hnRNP H3 is produced by E.coli expression system and the target gene encoding M1-Y346 is expressed with a His&SUMO tag at the N-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3189
UniProt ID:	P31942
Synonyms:	Heterogeneous nuclear ribonucleoprotein H3; Heterogeneous nuclear ribonucleoprotein 2H9; hnRNP 2H9; HNRPH3; HNRNPH3



[View online »](#)