

Product datasheet for **TP725101**

HNRPM (HNRNPM) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human hnRNP M
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Q158
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human hnRNP M is produced by E.coli expression system and the target gene encoding M1-Q158 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4670
UniProt ID:	P52272
Synonyms:	Heterogeneous nuclear ribonucleoprotein M; HNRPM
Protein Families:	Druggable Genome
Protein Pathways:	Spliceosome



[View online »](#)