

## Product datasheet for **TP725102**

### SF2 (SRSF1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SRSF1 Protein(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-T248
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SRSF1 Protein is produced by E.coli expression system and the target gene encoding M1-T248 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	6426
UniProt ID:	<a href="#">Q07955</a>



[View online »](#)