

## Product datasheet for **TP725103**

### **HNRNPA0 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human hnRNP A0(N-His-SUMO)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	M1-F305
<b>Tag:</b>	N-His-SUMO
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human hnRNP A0 is produced by E.coli expression system and the target gene encoding M1-F305 is expressed with a His&SUMO tag at the N-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	10949
<b>UniProt ID:</b>	<a href="#">Q13151</a>



[View online »](#)