

## Product datasheet for **TP725104**

### RCN2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Reticulocalbin-2/RCN2
Species:	Human
Expression cDNA Clone or AA Sequence:	G23-E209
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Reticulocalbin-2/RCN2 is produced by E.coli expression system and the target gene encoding G23-E209 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5955
UniProt ID:	<a href="#">Q14257</a>
Synonyms:	Calcium-binding protein ERC-55; E6-binding protein; E6BP; ERC55



[View online »](#)