

## Product datasheet for **TP725105**

### MVP Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Major vault protein/MVP
Species:	Human
Expression cDNA Clone or AA Sequence:	N544-R893
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Major vault protein/MVP is produced by E.coli expression system and the target gene encoding N544-R893 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9961
UniProt ID:	<a href="#">Q14764</a>
Synonyms:	Lung resistance-related protein; LRP
Protein Families:	Druggable Genome



[View online »](#)