

Product datasheet for **TP725108**

RHEB Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GTP-binding protein Rheb
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-M184
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GTP-binding protein Rheb is produced by E.coli expression system and the target gene encoding M1-M184 is expressed.
Stability:	12 months from date of despatch
Locus ID:	6009
UniProt ID:	Q15382



[View online »](#)