

## Product datasheet for **TP725108**

### **RHEB Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant Human GTP-binding protein Rheb
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-M184
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human GTP-binding protein Rheb is produced by E.coli expression system and the target gene encoding M1-M184 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6009
UniProt ID:	<u><a href="#">Q15382</a></u>



[View online »](#)