

Product datasheet for **TP725109**

PPP1R7 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PPP1R7
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-F360
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PPP1R7 is produced by E.coli expression system and the target gene encoding M1-F360 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5510
UniProt ID:	Q15435
Synonyms:	Protein phosphatase 1 regulatory subunit 7; Protein phosphatase 1 regulatory subunit 22
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)