

## Product datasheet for **TP725113**

### ASXL1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Polycomb group protein ASXL1
Species:	Human
Expression cDNA Clone or AA Sequence:	H61-R246
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Polycomb group protein ASXL1 is produced by E.coli expression system and the target gene encoding H61-R246 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	171023
UniProt ID:	<a href="#">Q8IXJ9</a>



[View online »](#)