

## Product datasheet for **TP725114**

### HDDC3 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human HDDC3/MESH1(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-E136
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human HDDC3/MESH1 is produced by E.coli expression system and the target gene encoding M1-E136 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	374659
UniProt ID:	<a href="#">Q8N4P3</a>



[View online »](#)