

## Product datasheet for **TP725118**

### **HNRPUL1 (HNRNPUL1) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human HNRNPUL1
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	M1-D227
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human HNRNPUL1 is produced by E.coli expression system and the target gene encoding M1-D227 is expressed.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	11100
<b>UniProt ID:</b>	<a href="#">Q9BUJ2</a>
<b>Synonyms:</b>	Heterogeneous nuclear ribonucleoprotein U-like protein 1; Adenovirus early region 1B-associated protein 5; E1B-55 kDa-associated protein 5; E1B-AP5; HNRPUL1
<b>Protein Families:</b>	Druggable Genome



[View online »](#)