

Product datasheet for **TP725119**

GPRC5C Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GPRC5C Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	S301-V433
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GPRC5C Protein is produced by E.coli expression system and the target gene encoding S301-V433 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55890
UniProt ID:	Q9NQ84
Synonyms:	G-protein coupled receptor family C group 5 member C; Retinoic acid-induced gene 3 protein; RAIG3
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)