

Product datasheet for **TP725120**

ATP6V1H Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human V-ATPase subunit H(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	S296-S483
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human V-ATPase subunit H is produced by E.coli expression system and the target gene encoding S296-S483 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	51606
UniProt ID:	Q9UI12



[View online »](#)