

Product datasheet for TP725125

alpha Adaptin (AP2A2) Human Recombinant Protein

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	Recombinant Human AP-2 complex subunit alpha-2/AP2A2
Species:	Human
Expression cDNA Clone or AA Sequence:	S590-F939
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human AP-2 complex subunit alpha-2/AP2A2 is produced by E.coli expression system and the target gene encoding S590-F939 is expressed.
Stability:	12 months from date of despatch
Locus ID:	161
UniProt ID:	<u>094973</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US