

Product datasheet for **TP725126**

PSMB3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Proteasome subunit beta type-3(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-D205
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Proteasome subunit beta type-3 is produced by E.coli expression system and the target gene encoding M1-D205 is expressed with a His&SUMO tag at the N-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5691
UniProt ID:	P49720
Synonyms:	Proteasome chain 13; Proteasome component C10-II; Proteasome theta chain; PSMB3
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Proteasome



[View online »](#)