

Product datasheet for **TP725127**

RPL11 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RPL11
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K178
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RPL11 is produced by E.coli expression system and the target gene encoding M1-K178 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6135
UniProt ID:	P62913
Synonyms:	60S ribosomal protein L11; CLL-associated antigen KW-12; Large ribosomal subunit protein uL5
Protein Pathways:	Ribosome



[View online »](#)