

Product datasheet for **TP725128**

GAB1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GRB2-associated-binding protein 1
Species:	Human
Expression cDNA Clone or AA Sequence:	T120-A400
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GRB2-associated-binding protein 1 is produced by E.coli expression system and the target gene encoding T120-A400 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2549
UniProt ID:	Q13480
Synonyms:	GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1; GAB1
Protein Families:	Druggable Genome
Protein Pathways:	ErbB signaling pathway, Neurotrophin signaling pathway, Renal cell carcinoma



[View online »](#)