

Product datasheet for **TP725134**

ARL8B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ARL8B(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	E20-S186
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ARL8B is produced by E.coli expression system and the target gene encoding E20-S186 is expressed with a His&SUMO tag at the N-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55207
UniProt ID:	<u>Q9NVJ2</u>



[View online »](#)