

Product datasheet for **TP725135**

Apc10 (ANAPC10) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human APC10
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-R185
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human APC10 is produced by E.coli expression system and the target gene encoding M1-R185 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10393
UniProt ID:	Q9UM13



[View online »](#)