

Product datasheet for **TP725136**

ANKRD42 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ANKRD42
Species:	Human
Expression cDNA Clone or AA Sequence:	Y225-N389
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ANKRD42 is produced by E.coli expression system and the target gene encoding Y225-N389 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	338699
UniProt ID:	Q8N9B4
Synonyms:	Ankyrin repeat domain-containing protein 42



[View online »](#)