

## Product datasheet for **TP725139**

### AUTS2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Autism susceptibility gene 2 protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A911-R1259
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Autism susceptibility gene 2 protein is produced by E.coli expression system and the target gene encoding A911-R1259 is expressed.
Stability:	12 months from date of despatch
Locus ID:	26053
UniProt ID:	<a href="#">Q8WXX7</a>



[View online »](#)