

Product datasheet for **TP725140**

ARFGEF2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Brefeldin A-inhibited GEP 2
Species:	Human
Expression cDNA Clone or AA Sequence:	I221-L370
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Brefeldin A-inhibited GEP 2 is produced by E.coli expression system and the target gene encoding I221-L370 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10564
UniProt ID:	Q9Y6D5
Synonyms:	Brefeldin A-inhibited guanine nucleotide-exchange protein 2; ADP-ribosylation factor guanine nucleotide-exchange factor 2; ARFGEF2; ARFGEP2; BIG2



[View online »](#)