

## Product datasheet for **TP725150**

### DIS3L2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human hDIS3L2
Species:	Human
Expression cDNA Clone or AA Sequence:	D27-F250
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human hDIS3L2 is produced by E.coli expression system and the target gene encoding D27-F250 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	129563
UniProt ID:	<a href="#">Q8IYB7</a>
Synonyms:	DIS3-like exonuclease 2; DIS3L2



[View online »](#)