

## Product datasheet for **TP725153**

### **DNAH5 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human DNAH5 Protein(N-His-SUMO)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	M1-D142
<b>Tag:</b>	N-His-SUMO
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human DNAH5 Protein is produced by E.coli expression system and the target gene encoding M1-D142 is expressed with a His&SUMO tag at the N-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	1767
<b>UniProt ID:</b>	<u><a href="#">Q8TE73</a></u>



[View online »](#)