

Product datasheet for **TP725157**

ELP3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human hELP3
Species:	Human
Expression cDNA Clone or AA Sequence:	W359-K547
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human hELP3 is produced by E.coli expression system and the target gene encoding W359-K547 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55140
UniProt ID:	Q9H9T3
Synonyms:	Elongator complex protein 3; tRNA uridine(34) acetyltransferase; ELP3



[View online »](#)