

Product datasheet for **TP725162**

FHDC1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FHDC1
Species:	Human
Expression cDNA Clone or AA Sequence:	R959-R1087
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human FHDC1 is produced by E.coli expression system and the target gene encoding R959-R1087 is expressed.
Stability:	12 months from date of despatch
Locus ID:	85462
UniProt ID:	<u>Q9C0D6</u>



[View online »](#)