

## Product datasheet for **TP725163**

### ACBD3 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Golgi resident protein GCP60
Species:	Human
Expression cDNA Clone or AA Sequence:	R229-R528
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Golgi resident protein GCP60 is produced by E.coli expression system and the target gene encoding R229-R528 is expressed.
Stability:	12 months from date of despatch
Locus ID:	64746
UniProt ID:	<u><a href="#">Q9H3P7</a></u>



[View online »](#)