

Product datasheet for **TP725165**

MOV10 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Helicase MOV-10
Species:	Human
Expression cDNA Clone or AA Sequence:	Q869-L1003
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Helicase MOV-10 is produced by E.coli expression system and the target gene encoding Q869-L1003 is expressed.
Stability:	12 months from date of despatch
Locus ID:	4343
UniProt ID:	Q9HCE1



[View online »](#)