

Product datasheet for **TP725166**

HOMER3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Homer protein homolog 3
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-P361
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Homer protein homolog 3 is produced by E.coli expression system and the target gene encoding M1-P361 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9454
UniProt ID:	Q9NSC5
Synonyms:	Homer-3; HOMER3
Protein Families:	Druggable Genome



[View online »](#)