

Product datasheet for **TP725171**

ISCA1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ISCA1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	S14-I129
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ISCA1 Protein is produced by E.coli expression system and the target gene encoding S14-I129 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	81689
UniProt ID:	Q9BUE6
Synonyms:	Iron-sulfur cluster assembly 1 homolog, mitochondrial; HESB-like domain-containing protein 2; Iron-sulfur assembly protein IscA; hlscA; HBLD2



[View online »](#)