

## Product datasheet for **TP725172**

### Nucleotide binding protein like (NUBPL) Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human Iron-sulfur protein NUBPL  |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | G39-E319   |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human Iron-sulfur protein NUBPL is produced by E.coli expression system and the target gene encoding G39-E319 is expressed.  |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 80224  |
| UniProt ID:                           | <a href="#">Q8TB37</a>   |
| Synonyms:                             | IND1 homolog; Nucleotide-binding protein-like; huInd1  |



[View online »](#)