

## Product datasheet for **TP725173**

### LMO7 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LMO-7
Species:	Human
Expression cDNA Clone or AA Sequence:	I696-Y838
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human LMO-7 is produced by E.coli expression system and the target gene encoding I696-Y838 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4008
UniProt ID:	<a href="#">Q8WWI1</a>
Synonyms:	LIM domain only protein 7; F-box only protein 20; LOMP; FBX20; FBXO20
Protein Pathways:	Adherens junction



[View online »](#)