

Product datasheet for **TP725174**

MTR Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Methionine synthase/MTR
Species:	Human
Expression cDNA Clone or AA Sequence:	R917-D1265
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Methionine synthase/MTR is produced by E.coli expression system and the target gene encoding R917-D1265 is expressed.
Stability:	12 months from date of despatch
Locus ID:	4548
UniProt ID:	<u>Q99707</u>



[View online »](#)