

Product datasheet for **TP725176**

MTUS1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MTUS1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1031-R1270
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human MTUS1 Protein is produced by E.coli expression system and the target gene encoding M1031-R1270 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	57509
UniProt ID:	Q9ULD2
Synonyms:	Microtubule-associated tumor suppressor 1; AT2 receptor-binding protein; Angiotensin-II type 2 receptor-interacting protein; Mitochondrial tumor suppressor 1;ATBP; ATIP; GK1; MTSG1
Protein Families:	Druggable Genome, Transcription Factors, Transmembrane



[View online »](#)