

## Product datasheet for **TP725178**

### MDFI Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MDFI Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Q160
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human MDFI Protein is produced by E.coli expression system and the target gene encoding M1-Q160 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4188
UniProt ID:	<a href="#">Q99750</a>
Synonyms:	MyoD family inhibitor; Myogenic repressor I-mf
Protein Families:	Transcription Factors



[View online »](#)