

Product datasheet for **TP725179**

MYO3A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myosin-IIIa/ MYO3A
Species:	Human
Expression cDNA Clone or AA Sequence:	Q1090-T1307
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Myosin-IIIa/ MYO3A is produced by E.coli expression system and the target gene encoding Q1090-T1307 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	53904
UniProt ID:	Q8NEV4
Synonyms:	None
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)