

Product datasheet for **TP725181**

Neurobeachin (NBEA) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Neurobeachin/NBEA
Species:	Human
Expression cDNA Clone or AA Sequence:	T1589-L1703
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Neurobeachin/NBEA is produced by E.coli expression system and the target gene encoding T1589-L1703 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	26960
UniProt ID:	Q8NFP9
Synonyms:	Lysosomal-trafficking regulator 2; Protein BCL8B; BCL8B; LYST2



[View online »](#)