

Product datasheet for **TP725187**

PTRH2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PTH 2
Species:	Human
Expression cDNA Clone or AA Sequence:	G63-Y179
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PTH 2 is produced by E.coli expression system and the target gene encoding G63-Y179 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51651
UniProt ID:	Q9Y3E5
Synonyms:	Peptidyl-tRNA hydrolase 2, mitochondrial; Bcl-2 inhibitor of transcription 1; BIT1; PTRH2
Protein Families:	Transmembrane



[View online »](#)