

## Product datasheet for **TP725190**

### PRAM1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PRAM-1
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-E299
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human PRAM-1 is produced by E.coli expression system and the target gene encoding M1-E299 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	84106
UniProt ID:	<a href="#">Q96QH2</a>
Synonyms:	PML-RARA-regulated adapter molecule 1; PRAM
Protein Families:	Druggable Genome



[View online »](#)