

Product datasheet for **TP725192**

Prickle (PRICKLE1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Prickle-like protein 1
Species:	Human
Expression cDNA Clone or AA Sequence:	A321-M580
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Prickle-like protein 1 is produced by E.coli expression system and the target gene encoding A321-M580 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	144165
UniProt ID:	Q96MT3
Synonyms:	REST/NRSF-interacting LIM domain protein 1; PRICKLE1; RILP
Protein Families:	Druggable Genome
Protein Pathways:	Wnt signaling pathway



[View online »](#)