

## Product datasheet for **TP725193**

### **BICC1 Human Recombinant Protein**

#### **Product data:**

|  |  |
|--|--|
| <b>Product Type:</b>                         | Recombinant Proteins   |
| <b>Description:</b>                          | Recombinant Human Protein bicaudal C homolog 1   |
| <b>Species:</b>                              | Human  |
| <b>Expression cDNA Clone or AA Sequence:</b> | M1-S130  |
| <b>Buffer:</b>                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| <b>Note:</b>                                 | MADE TO ORDER   Recombinant Human Protein bicaudal C homolog 1 is produced by E.coli expression system and the target gene encoding M1-S130 is expressed.  |
| <b>Storage:</b>                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| <b>Stability:</b>                            | 12 months from date of despatch  |
| <b>Locus ID:</b>                             | 80114  |
| <b>UniProt ID:</b>                           | <a href="#">Q9H694</a>   |
| <b>Synonyms:</b>                             | BICC1; Bic-C   |



[View online »](#)