

## Product datasheet for **TP725197**

### **BOULE (BOLL) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Protein boule-like
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	M1-Y283
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Protein boule-like is produced by E.coli expression system and the target gene encoding M1-Y283 is expressed.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	66037
<b>UniProt ID:</b>	<a href="#">Q8N9W6</a>
<b>Synonyms:</b>	BOLL; BOULE



[View online »](#)