

Product datasheet for **TP725198**

PGEA1 (CBY1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein chibby homolog 1/CBY1(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K126
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein chibby homolog 1/CBY1 is produced by E.coli expression system and the target gene encoding M1-K126 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	25776
UniProt ID:	Q9Y3M2



[View online »](#)