

Product datasheet for **TP725200**

Slingshot homolog 1 (SSH1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase Slingshot homolog 1/SSH1
Species:	Human
Expression cDNA Clone or AA Sequence:	M928-S1049
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein phosphatase Slingshot homolog 1/SSH1 is produced by E.coli expression system and the target gene encoding M928-S1049 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54434
UniProt ID:	Q8WYL5
Synonyms:	SSH-like protein 1; SSH-1L; hSSH-1L
Protein Families:	Druggable Genome, Phosphatase
Protein Pathways:	Regulation of actin cytoskeleton



[View online »](#)