

Product datasheet for **TP725202**

LZTR2 (SEC16B) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein transport protein Sec16B
Species:	Human
Expression cDNA Clone or AA Sequence:	Q48-M273
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein transport protein Sec16B is produced by E.coli expression system and the target gene encoding Q48-M273 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	89866
UniProt ID:	Q96JE7
Synonyms:	Leucine zipper transcription regulator 2; Regucalcin gene promoter region-related protein p117 (RGPR-p117); SEC16 homolog B



[View online »](#)